

CADD SOLIS INFUSION PUMP INSTRUCTIONS

1. Wash hands and gather supplies.
 - a. Alcohol Prep Pads
 - b. 2 Saline syringes
 - c. Swab cap
 - d. 1 Heparin syringe
 - e. Pump
 - f. Tubing
 - g. Solution (TPN/Fluids)
2. Power on the pump.
3. Pull the latch open on the pump, attach tubing and close latch as demonstrated by your healthcare professional.
4. When the solution (medication/TPN/fluids) is ready, spike the container.
5. Follow the prompt on the pump, select yes to prime (MAKE SURE PATIENT IS NOT CONNECTED TO TUBING AT THIS TIME).
6. When there is a consistent drip coming out of the tubing, select "Stop Priming" on the pump and select "No" on continue priming (unless priming is still needed).
7. Select "Review" on pump and review all settings. If pump is prompting a code please use #201. Compare settings to the label on the solution, if correct select "Accept Value" to each step.
8. Remove Swab Cap from PICC line paying close attention not to contaminate the clave. If at any time you contaminate the end of your line, clean it with an alcohol prep pad with friction for 15 seconds.
9. Attach Saline syringe and flush using a pulsing motion. NOTE - if there is resistance -STOP- and call your Home Health Nurse.
10. Attach tubing to patient and unclasp clamp.
11. Push "Yes" on pump to start.
12. It is important to "Reset Reservoir Volume" when hooking up to a new bag.
13. When infusion is complete, disconnect solution from PICC line.
14. Attach Saline syringe and flush with a pulsing motion (do not flush or use Heparin when Milrinone or Dobutamine is running continuous).
15. Attach Heparin syringe and flush with a pulsing motion.
16. Cap off with a new Swab Cap.

You can access the expertise of a pharmacist round the clock, every day of the week. Should you have any inquiries, concerns, or require consultation, please don't hesitate to reach out at **(530) 225-8898**. Our team is available at all times to assist you.

For your safety, your medication is on a schedule as prescribed by your clinician. To ensure proper delivery, we will always contact you prior to sending your next dose. It is your responsibility to contact us when you are low on supplies. We look forward to hearing from you to request more saline, heparin, PICC kits, alcohol prep pad and / or swab caps. Your safety is our greatest concern. Wellscrip (530) 225-8898.