

24 HOURS INFUSOR INSTRUCTIONS

1. Wash hands
2. Gather supplies
 - a. Alcohol prep pads
 - b. 1 Saline syringe
 - c. Infusor
3. Clean the end of the clave with an alcohol prep pad with friction for 15 seconds. When working with your PICC line pay close attention not to contaminate the clave. If at any time you contaminate the end of your line repeat cleaning step.
4. Attach saline syringe and flush using a pulsing method.
Note: if there is resistance STOP and call your Home Health Nurse.
5. Attach infusor and unclasp clamp.
6. When the infusion is completed 24 hours later, detach the infusor.
7. Start over again at Step 1.

For your safety, your medication is on a schedule as prescribed by your clinician. To ensure proper delivery, we will always contact you prior to sending your next dose. It is your responsibility to contact us when you are low on supplies. We look forward to hearing from you to request more saline, heparin, PICC kits, alcohol prep pad and/or swab caps. Your safety is our greatest concern.

You can access the expertise of a pharmacist 24/7. Should you have any inquiries, concerns, or require consultation, contact us at **(530) 225-8898**. Our team is available at all times to assist you.

Thank you.